

# A GAME-CHANGING RNA DELIVERY PLATFORM COMBINING HIGH CELL ENTRY AND MULTIPLE RNA TRANSFER FOR NEXT GENERATION THERAPY : LENTIFLASH®

flash  
THERAPEUTICS

FLORINE SAMAIN<sup>1</sup>, ALEXANDRA ICHÉ<sup>1</sup>, NICOLAS MARTIN<sup>1</sup>, CHRISTINE DUTHOIT<sup>1</sup>, RÉGIS GAYON<sup>1</sup>, EMELINE BENUZZI<sup>2</sup>, BARBARA GARMY-SUSIN<sup>2</sup> AND PASCALE BOUILLÉ<sup>1</sup>


<sup>1</sup>FLASH THERAPEUTICS, TOULOUSE, FRANCE, <sup>2</sup>INSTITUTE OF METABOLIC AND CARDIOVASCULAR DISEASES OF TOULOUSE (I2MC), INSERM UMR 1297, CHU RANGUEIL, 1 AV JEAN POULHÈS, 31432 TOULOUSE CEDEX 4, FRANCE

## A. Challenges of RNA therapies

**Gene therapy** approaches show that there is no universal delivery tool for all therapeutic strategies. Compared to DNA delivered-therapies, **RNA therapies** are expected to be more versatile, cover a broad range of applications with minimal regulatory concerns and thus address a large variety of diseases. The technology targets applications in which a **transient expression** is expected.

As a game-changing RNA carrier, **LentiFlash®**, a non-integrative bacteriophage-lentivirus chimera, can efficiently and safely deliver **multiple RNA species** that are **rapidly bioavailable**, leading to a **high and short-term expression** of the transferred messenger into the cell cytoplasm.

### LENTIFLASH® PARTICLE



Reference : (Prel et al. Mol Ther Methods Clin Dev. 2015)

## C. Clinical Context and mouse model

**Lymphedema** is a lymphatic vascular system disorder characterized by:

- impaired lymphatic return and swelling of the extremities
- accumulation of undrained interstitial fluid/lymph.


It results in fibrosis and adipose tissue deposition in the affected area. It can occur after cancer surgery and lymph node removal. Indeed, 10-15% of women develop arm secondary lymphedema after surviving breast cancer. There is no curative treatment for lymphedema.

**TherAlymph** clinical program's main objective is to establish a multiple gene therapy for secondary lymphedema, based on the transient expression of two factors allowing the restoration of a normal lymphatic function. LentiFlash® will be used to deliver the two therapeutic RNAs by intra-dermal injections.




### LYMPHEDEMA MOUSE MODEL

Secondary lymphedema is induced in mice by mastectomy of the 2<sup>nd</sup> mammary gland, associated with brachial and axillary lymphadenectomy. LentiFlash® (LF) is intra-dermally injected in limb at day 10.



## B. What is LentiFlash® ?




This **biological RNA delivery technology** mediated by a lentiviral particle is an attractive approach as it combines most of the inherited properties of lentiviral vectors (cell entry and tropism) without the potential adverse effects from long-lasting expression or genomic integration.

The lentiviral external structure (membrane, capsid) ensures **RNAs protection**.

➔ **Scalable production process already validated for clinical settings**

### RNA DELIVERY FOR A TRANSIENT EXPRESSION



- No GMO generated
  - No reverse transcription
  - No DNA integration, no insertional mutagenesis
- Rapid & transient expression

## D. Pre-clinical results

After lymph nodes resection, mice exhibit a reproducible reduction of lymphatic drainage associated with dermal backflow (fig. 1) and increase of leg diameter (fig. 2) up to 4 weeks post-surgery.




Fig. 1: Lymphographies of treated or untreated limbs from control and lymphedema mice at day 14 post-surgery




Fig. 2: Measure of leg diameter

### Lymphedema resolution after therapeutic RNAs delivery by LentiFlash®

- Lymphatic vessels restoration and suppression of dermal backflow (fig. 1)
- Decrease of limb diameter upon LF injection at day 10, reaching the normal control value at day 21 (fig. 2)

### LentiFlash® induces a strong protein expression increase

At day 14 after injection of LF expressing Apelin+VEGF-C or Apelin alone, or an integrative Lentivector (iLV) expressing both genes, we observed a 2-fold increase of Apelin expression level, showing that LF is able to provide equivalent protein expression as iLV *in vivo*, even when it delivers 2 different mRNAs (fig. 3).




Fig. 3: Quantification of Apelin expression in plasma

## E. Conclusion

The LentiFlash® properties, associated with our own **lentiviral production platform** compliant with the **cGMPs**, offer additional safety considerations making it a versatile and safe mean for human therapy.

A **first-in-human phase I/IIa clinical trial** on patients who developed lymphedema after breast cancer using RNA delivery by LentiFlash® will be performed in 2024 at the Toulouse University Hospital, France.

Furthermore, LentiFlash®, as an RNA delivery tool, can be used for a broad range of applications, such as **gene editing** (Mianné et al. 2022) and **vaccination/immunotherapy applications** for both infectiology and oncology purposes.

- LentiFlash® packages **biological RNAs of human origin**, without any lentiviral sequence.
- LentiFlash® can deliver **multiple RNAs**.
- LentiFlash® does **not** lead to **adverse immune responses**.
- Transduction by LentiFlash® does **not** result in **GMO generation**.
- It combines the efficient delivery of lentiviral vectors with the safety of RNA delivery since it enables **highly efficient transfer** and **transient expression**:

- ✓ LentiFlash® particles display a very large tropism thanks to **VSV-G pseudotyping**.
- ✓ **All cell types** can be efficiently transduced by LentiFlash® without altering cell viability nor phenotype.

